

Pressure Mat Switch

1: Switch Type

SIZE

MAC-720	540mm*780mm
MAC-721	300mm*790mm
MAC-722	300mm*300mm
MAC-160	230mm*160mm
MAC-170	595mm*170mm
MAC-305	305ø*4mm
MAC-350	350mm*500mm
MAC-390	720mm*390mm
MAC-540	540mm*395mm
MAC-560	720mm*560mm
MAC-720RB (Yoga)	550mm*800mm
MAC-1200RB (Yoga)	600mm*1200mm
MAC-1600RB (Yoga)	600mm*1600mm

2: PVC Color + Use

B (Black)	Floor Use (1 Year Guarantee)
WF (White)	Floor Use (1 Year Guarantee)
WM6M (White+Medical)	Medical Use (6 Month Guarantee)
WM12M (White+Medical)	Medical Use (1 Year Guarantee)

Guiding of model:

MAC-720 - **B** - **LMP** - **150** - **2 Wires**

1 2 3 4 5

3: Leadwire Type

L	Leadwire Only
MP	6.3 male mono plug
LMP	6.3 L type male mono plug
FMP	6.3 female mono plug
IMP	3.5 male mono plug
ILMP	3.5 L type male mono plug
IFMP	3.5 female mono plug
JMP	2.5 male mono plug
JLMP	2.5 L type male mono plug
JFMP	2.5 female mono plug
RJ-11	TEL PLUG

4: Leadwire Length (cm)

30	100	200
70	150	500

5: Loop Type

2 Wires
4 Wires



MULTICE & Co., LTD.
彰信企業股份有限公司

P.O. BOX 47-91 TAIPEI, TAIWAN. 111 台北市士林區中正路668號3F

TEL: 886-2-28123313 FAX: 886-2-28123316

http://www.multice.com e-mail: multice@ms34.hinet.net

Visitor chime (MVT)

Place the pressure mat sensor underneath a larger size floor mat in front of an entrance, at the top or bottom of stairs, in front of a cabinet or room, next to a bed, anywhere that you need to be notified if someone is moving or has arrived. When someone steps on the pressure mat, the alarm unit will sound for approximately 30 seconds.



MVT Alarm is ideal for an entrance farther apart. Features include:

- Always "ON" for resident safety
- Flashing alert light will instantly identify roommate in need...day or night
- Volume choice of 90 or 75 dB
- Three AAA batteries included
- Low battery indicator light
- Test & reset button
- AC/DC adapter option (sold separately)
- Alarm warning can be delayed by 30sec
- Can be used with pressure mat (EX: MAC-720-B-JLMP-150...)

Guiding of model: **MVT** - **1000** - **W** - **NO** - **1**

1: Model

MVT	Visitor chime
MFT	Fall Prevention Alarm

2: Product number

1000	cord
2000	wireless

3: Color

W	White
B	Black

4: Input signal

NO
NC

5: Output signal

1	None
2	OC output

Fall Prevention Alarms (MFT)

This alarm is to keep an eye on the patient who is on a wheelchair or a sick bed. Once the patient departs from the chair or bed pad, the alarm be activated for calling nurse to help.

The most simple fall prevention alarm that can be used for retirement communities as well as family homes.

MFT Alarm is ideal for rooms farther apart. Features include:

- Always "ON" for resident safety
- Flashing alert light will instantly identify roommate in need...day or night
- Volume choice of 90 or 75 dB
- Three AAA batteries included
- Low battery indicator light
- Test & reset button
- AC/DC adapter option (sold separately)
- Can be used with bed mat (EX: MAC-720-WM6M-JLMP-150...), chair pad (EX: MAC-722-WM6M-JLMP-150...), nurse cord (EX: EN-EH110-NO...)

